Rotating Anode X-Ray Tube Tubes
Radiogènes à Anode Tournante
Drehanoden - Röntgenröhre
Tubos de Rayos-X con Ánodo Giratorio

Note: Document originally drafted in the English language.

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

Le Gaine MX165NP

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Encaje de MX165NP

Maximum Peak Voltage	140 kV
Anode to Ground	70 kV
Cathode to Ground	70 kV
Maximum X-ray Tube Assembly Heat Content	2.6 MJ (3.5 MHU)
Maximum Continuous Heat Dissipation (max. housing temperature 78°C) (Includes stator heat)	3.7 kW (5.0 kHU/sec)
Focal Point Position (Central Ray) Within 1mm (X, Y Direction from the center of radiation port.)	
X-Ray Tube Assembly	
Permanent filtration	1.2 mm Al IEC 60522
Loading Factors for Leakage Radiation	140 kV, 26 mA
Ambient Air Temperature Limits for Operation	5°C to 35°C
Temperature Limits for Storage and Transport	-34°C to +60°C
Humidity	+10% to +90%
Atmospheric Pressure Range	70 kPa to 106 kPa
Weight - Housing	70 kg (154 lbs)
IEC Classification	Class 1
Safety Devices - Thermal Switch	
Normally Closed Contact	Opening at 74°C to 79°C
Pressure Switch	
Normally Closed Contact	Opening at 5 PSIG ±3 PSIG

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
GS-3576S TUBE HOUSING ASSEMBLY HEATING AND COOLING

Note:
 1. Heat inputs into housing include tube power, filament power, and stator power.
 2. Heating curves based on no restrictions of natural convection around tube housing assembly.
 3. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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

Housing Connection Stator and Safety Switches Plug Connector

Heat Exchanger Connection Pump and Fans 115 VAC, 60 Hz Plug Connector

Stator Drive Frequency Fréquence d'entraînement du stator Statorantrieb Frequenz Frecuencia de la impulsión del estator	RPM
50 Hz	2800 - 3000
60 Hz	3400 - 3600
150 Hz	8500 - 9000
180 Hz	9500 - 10,800


Stator:		
Stator Coil Resistance:		
Black to White	23Ω ±15%	
Green to White	23Ω ±15%	
Black to Green	46Ω ±15%	
Starter Voltage:		
	Start	Run
50/60 Hz	385 VAC	75 VAC
150/180 Hz	450 VAC	85 VAC
Time to Full Speed:		
50/60 Hz	0 - 2700 RPM 10 Sec.	
X-Ray Tube Assembly:		
GS-3576S/MX165NP IEC 60601-2-28		

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ANODE HEATING AND COOLING CURVES

Note:
 1. Heating and cooling curves reflect maximum tube performance. Tube operation is ultimately limited by system software control.

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