

# Neurovirtual EEG

Neurovirtual is the best solution for your Neurology Department, no matter what your needs are—Mobile Routine EEG or Full Long-Term Monitoring.

## Easy to Use



The Neurovirtual BWAnalysis is the most user friendly and intuitive EEG software on the market. Everything from setting up to performing an exam to reading EEG studies is easy, with all icons displayed right on the screen.

## Free Reading Station

The BWII can be used in an active network environment, enabling you to read your EEG results from any computer you choose. You can read from the computer originally used for the test, your office, or even the comfort of your own home, free of charge.

## Remote Access

Studies can be accessed by an authorized physician from any computer or Apple iPad\* to review the EEG data, even during the live recording\*\*.



## BWAnalysis EEG Highlights

- Brain Mapping Included on the Analysis Software
- Works with most EMR platforms and HL7 interface
- Custom reports in Microsoft Word for each physician and practice
- User-designed workspaces for collection, analysis and MD Review
- Montages, filters, sensitivity and page speed can be adjusted during recording and reading
- Plug-n-play interface for spike detection software

\*Apple iPad is a trademark of Apple Computer \*\*A 3rd party app will be used for the live review interface


## Neurovirtual Amplifiers

Specifications	BWII EEG	BW3 EEG
Total Channels	25 Channels	50 Channels
DC Channels	NA	8 DCX
Communication	Ethernet	Ethernet
Studies Performed	EEG Only	EEG and PSG
Photic Stimulator	White Led 1-30Hz	White Led 1-30Hz
Impedance Check	Software	Software and Bedside
Notch Filter	50/60Hz	50/60Hz
Low Filters	0.03Hz to 15Hz	0.03Hz to 15Hz
High Filters	1Hz to 100 Hz	1Hz to 100 Hz

### BWII EEG

#### *Lightweight - Portable - Versatile*

Portable equipment that can be transported securely in a suitcase to different locations or stored in a stationary EEG cart.



### BW3 EEG

#### *The Complete System*

- All inputs found on the headbox are duplicated on the amplifier, to offer EEG professionals options they never had before
- Bedside Impedance Check with LED indicators directly on the amplifier

Contact us today for a product demonstration:  
 +1-786-693-8200 or  
[info@neurovirtual.com](mailto:info@neurovirtual.com)